



Repose Range Overview



Reactive Pressure Redistribution

Simple yet clinically effective

Healthcare professionals involved in the treatment and prevention of pressure ulcers are often faced with time constraints, limited resources and increasing budget pressures all of which can impact on patient outcomes.

Repose® is a range of static air, pressure redistribution and reduction support surfaces, all with high levels of immersion and envelopment. A combination of unique materials and Smart Valve Technology, has been clinically proven to deliver statistically significant patient benefits.¹



Repose protects patients at all levels of risk



Unlike, other static air products and dynamic systems, through Randomized Control Trials, Repose reduces peak and average pressures and consistently delivers low levels of pressure ulcer incidence, in a cost effective, easy to use system.¹⁻⁵

Repose is indicated for the prevention of pressure ulcers for patients at very high risk and treatment of all categories of pressure ulcer when used as part of a package of care.¹⁻⁵



Repose is indicated for the treatment of all categories of pressure ulcers¹⁻⁵

Benefits of Repose



Clinically effective

Repose has demonstrated low levels of pressure ulcer incidence across multiple studies (5.2%¹, 4.8%², 5.1%³, 5.2%⁴)



Smart Valve Technology

Easy and reliable, Repose Smart Valve Technology ensures optimal inflation.



Durable

If cared for correctly, Repose will provide pressure ulcer prevention for several years⁵.



Support

Dedicated sales and clinical support ensure clinicians and patients are supported every step of the care pathway



Comfortable

Repose is breathable and has no sharp edges, increasing comfort, compliance and prevents tissue damage.



2 year warranty

Repose is supplied with a 2 year warranty as standard

The Repose Range



Surface Solutions



Repose® Mattress Overlay

Order code: 6011100

Dimensions: L 74in x W 30in
L 190cm x W 77cm

Weight Limit: 489lbs*



Repose® Mattress Overlay & Cushion Set

Order code: 6121101

Dimensions: Mattress L 74in x W 30in
L 190cm x W 77cm
Cushion L 18in x W 18in
L 45cm x W 45cm

Weight Limit: Cushion - 306lbs* Mattress - 489lbs*



Seating Solutions



Repose® Cushion

Order code: 6201100

Dimensions: L 18in x W 18in
L 45cm x W 45cm

Weight Limit: 306lbs



Repose® Care-Sit

Order code: 6351001

Dimensions: L 37in x W 18in
L 95cm x W 45cm

Weight Limit: 306lbs



Foot Care Solutions

New Sizes



Repose® Foot Protector Standard

Order code: 6501100

Dimensions: L 15in x W 8in x H 9in
L 38cm x W 20cm x H 23cm

Weight Limit: 306lbs

Available sizes: XS S STD L

New Sizes



Repose® Foot Protector Plus Standard

Order code: 6511100

Dimensions: L 15in x W 8in x H 9in
L 38cm x W 20cm x H 23cm

Weight Limit: 306lbs

Available sizes: XS S STD



Repose® Trolley Mattress Overlay

Order code: 6601100
Dimensions: L 70in x W 22in
L 178cm x W 55cm
Weight Limit: 306lbs



Repose® Companion

Order code: 6990003
Dimensions: L 70in x W 22in
L 197cm x W 88cm
Weight Limit: 306lbs



Repose® Contur Overlay

Order code: 6310001
Dimensions: L 69in x W 21in
L 176cm x W 52cm
Weight Limit: 306lbs



PU Strap

Repose® Wedge

Order code: 6711100
Dimensions: L 27in x W 18in
L 69cm x W 45cm
Weight Limit: 306lbs

Repose® Wedge With PU Strap

Order code: 6711111



new

Repose® Flex Standard

Order code: 6750002
Dimensions: L 20in x W 28in
L 50cm x W 70cm
Weight Limit: 306lbs

Available sizes: XS S STD



Repose® Trolley Companion

Order code: 6990050

Dimensions: L 79in x W 24in
L 200cm x W 61.5cm

Weight Limit: 306lbs



Paediatric Solutions



Repose® Paediatric Mattress Overlay

Order code: 6801100

Dimensions: L 51in x W 26in
L 129cm x W 67cm

Weight Limit: 306lbs



Repose® Babynest

Order code: 6801101

Dimensions: L 23in x W 13in x H 4in
L 58cm x W 33cm x H 10cm

Weight Limit: 7lbs



Repose® Sole Protector

Order code: 6930001

Dimensions: L 34in x H 15in
L 86cm x H 37cm

Weight Limit: 306lbs



Repose® Babytherm


Order code: 6401100

Dimensions: L 30in x W 21in
L 76cm x W 53cm

Weight Limit: 306lbs


Accessories



 **Repose® Cushion Cover**
Order code: 6211100

 **Repose® Mattress Cover**
Order code: 6121101



 **Repose® Cushion Strap Extender**
Order code: 6000001
Dimensions: L 39in x W 1in
L 100cm x W 2.5cm

References

1. Beeckman et al. (2018) A multi-centre prospective randomised controlled clinical trial to compare the effectiveness and cost of a static air mattress and alternating air pressure mattress to prevent pressure ulcers in high-risk nursing home residents. *International Journal of Nursing Studies*, June, <https://doi.org/10.1016/j.ijnurstu.2019.05.015>
2. Van Leen, M. Hovius, S. Neyens, J. Hafens, R. Shols, J. (2011) Pressure relief, cold foam or static air? A single centre, prospective, controlled randomized clinical trial in a Dutch nursing home. *Journal of Tissue Viability*, 20 (1), 30-34
3. Beeckman et al. (2016) Static air support surfaces to prevent pressure injuries. *J Wound Ostomy Continence Nurs.* 2016;43(4):375-378
4. Van Leen et al. (2013) Pressure Relief with Visco-Elastic Foam or with Combined Static Air Overlay? A prospective, Crossover Randomized Clinical Trial in a Dutch Nursing Home. *Wounds* 2013;25(10):287-292
5. Poole M. (2009) A long-term durability assessment of static air-filled pressure redistribution devices in the community setting. *Wounds UK*.

*Mattress overlay when used on a high specification foam mattress 489lbs.

*Mattress overlay when used on a non-high specification foam mattress 306lbs.

For more information:

Call (800) 303-9544

Email: contactus@frontier-therapeutics.com

Website: www.frontier-therapeutics.com



300 International Drive, Suite 100 Buffalo,
NY 14221 United States of America

 © Copyright Frontier Therapeutics Limited, 2019
All rights reserved. GP0041_A70_02

repose[®]